# Undersøgelse for Carbapenemase Producerende Organismer (CPO)

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## Referencer

1. Girlich D, Bouihat N, Poirel L, Benouda A, Nordmann P. High rate of faecal carriage of extended-spectrum beta-lactamase and OXA-48 carbapenemase-producing Enterobacteriaceae at a university hospital in Morocco. *Clin Microbiol Infect* 2014;20(4):350-4.
2. Papadimitriou-Olivgeris M, Bartzavali C, Christofidou M, Bereksi N, Hey J, Zambardi G, et al. Performance of chromID(R) CARBA medium for carbapenemase-producing Enterobacteriaceae detection during rectal screening. *Eur J Clin Microbiol Infect Dis* 2014;33(1):35-40
3. Girlich D, Poirel L, Nordmann P. Comparison of the SUPERCARBA, CHROMagar KPC, and Brilliance CRE screening media for detection of Enterobacteriaceae with reduced susceptibility to carbapenems. *Diagn Microbiol Infect Dis* 2013;75(2):214-7.
4. Girlich D, Anglade C, Zambardi G, Nordmann P. Comparative evaluation of a novel chromogenic medium (chromID OXA-48) for detection of OXA-48 producing Enterobacteriaceae. *Diagn Microbiol Infect Dis* 2013;77(4):296-300.
5. Wilkinson KM, Winstanley TG, Lanyon C, Cummings SP, Raza MW, Perry JD. Comparison of four chromogenic culture media for carbapenemase-producing Enterobacteriaceae. *J Clin Microbiol* 2012;50(9):3102-4.
6. Vrioni G, Daniil I, Voulgari E, Ranellou K, Koumaki V, Ghirardi S, et al. Comparative evaluation of a prototype chromogenic medium (ChromID CARBA) for detecting carbapenemase-producing Enterobacteriaceae in surveillance rectal swabs. *J Clin Microbiol* 2012 Jun;(6):1841-6.
7. Singh K, Mangold KA, Wyant K, Schora DM, Voss B, Kaul KL, Hayden MK, Chundi V, Peterson LR. Rectal screening for *Klebsiella pneumoniae* carbapenemases: comparison of real-time PCR and culture using two selective screening agar plates. *J Clin Microbiol*. 2012;50(8):2596-600.
8. Adler A, Navon-Venezia S, Moran-Gilad J, Marcos E, Schwartz D, and Carmeli Y. Laboratory and Clinical Evaluation of Screening Agar Plates for Detection of Carbapenem-Resistant *Enterobacteriaceae* from Surveillance Rectal Swabs. *J Clin Microbiol*. 2011; 49(6): 2239–2242.
9. Samra Z, Bahar J, Madar-Shapiro L, Aziz N, Israel S, and Bishara J. Evaluation of CHROMagar KPC for Rapid Detection of Carbapenem-Resistant *Enterobacteriaceae*. *J Clin Microbiol*. 2008; 46(9): 3110–3111.
10. Diederer B, Chang C, Euser S, Stuart JC. Evaluation of four screening protocols for detection of extended-spectrum beta-lactamase-producing members of the Enterobacteriaceae. *J Med Microbiol* 2012;61(3):452-3.
11. Murk JL, Heddema ER, Hess DL, Bogaards JA, Vandenbroucke-Grauls CM, Debets-Ossenkopp YJ. Enrichment broth improved detection of extended-spectrum-beta-lactamase-producing bacteria in throat and rectal surveillance cultures of samples from patients in intensive care units. *J Clin Microbiol* 2009;47(6):1885-7.
12. Hansen F, Hammerum AM, Skov R, Haldorsen B, Sundsfjord A, Samuelsen O. Evaluation of the total MBL confirm kit (ROSCO) for detection of metallo- $\beta$ -lactamases in *Pseudomonas aeruginosa* and *Acinetobacter baumannii*. *Diagn Microbiol Infect Dis*. 2014;79(4):486-8.
13. Girlich D, Poirel L, Nordmann. Value of the Modified Hodge Test for Detection of Emerging Carbapenemases in Enterobacteriaceae. *J Clin Microbiol*. 2012; 50(2): 477–479.
14. Wang P, Chen S, Guo Y, Xiong Z, Hu F, Zhu D, Zhang Y. Occurrence of False Positive Results for the Detection of Carbapenemases in Carbapenemase-Negative *Escherichia coli* and *Klebsiella pneumoniae* Isolates. *PLoS One*. 2011;6(10).
15. Carvalhaes CG, Picão RC, Nicoletti AG, Xavier DE, Gales AC. Cloverleaf test (modified Hodge test) for detecting carbapenemase production in *Klebsiella pneumoniae*: be aware of false positive results *J Antimicrob Chemother*. 2010;65(2):249-51

16. Pasteran FG, Mendez T, Guerriero L et al. Sensitive screening tests for suspected class A carbapenemase production in species of Enterobacteriaceae. *J Clin Microbiol* 2009; 47: 1631–9.
17. Nordmann P, Poirel L, Dortet L Rapid detection of carbapenemase-producing Enterobacteriaceae. *Emerg Infect Dis.* 2012;18(9):1503-7.
18. Dortet L, Poirel L, Nordmann P. Rapid identification of carbapenemase types in Enterobacteriaceae and *Pseudomonas* spp. by using a biochemical test. *Antimicrob Agents Chemother.* 2012;56 (12):6437-40.
19. Vasoo S, Cunningham SA, Kohner PC, Simner PJ, Mandrekar JN, Lolans K, Hayden MK, Patel R. Comparison of a novel, rapid chromogenic biochemical assay, the Carba NP test, with the modified Hodge test for detection of carbapenemase-producing Gram-negative bacilli. *J Clin Microbiol.* 2013;51(9):3097-101.
20. Milillo M, Kwak YI, Snesrud E, Waterman PE, Lesho E, McGann P. Rapid and simultaneous detection of blaKPC and blaNDM by use of multiplex real-time PCR. *J Clin Microbiol.* 2013;51(4):1247-9.
21. Dortet L, Poirel L, Nordmann P. Further proofs of concept for the Carba NP test. *Antimicrob Agents Chemother.* 2014;58(2):1269
22. Hrabák J, Walková R, Studentová V, Chudácková E, Bergerová T. Carbapenemase activity detection by matrix-assisted laser desorption ionization-time of flight mass spectrometry. *J Clin Microbiol.* 2011;49(9):3222-7.
23. Burckhardt I, Zimmermann S. Using Matrix-Assisted Laser Desorption Ionization-Time of Flight Mass Spectrometry To Detect Carbapenem Resistance within 1 to 2.5 Hours. *J Clin Microbiol.* 2011; 49(9): 3321–3324.
24. Hrabák J, Studentová V, Walková R, Zemlicková H, Jakubu V, Chudácková E, Gniadkowski M, Pfeifer Y, Perry JD, Wilkinson K, Bergerová T. Detection of NDM-1, VIM-1, KPC, OXA-48, and OXA-162 carbapenemases by matrix-assisted laser desorption ionization-time of flight mass spectrometry. *J Clin Microbiol* 2012;50(7):2441-3.
25. Tzouvelekis LS, Markogiannakis A, Psychogiou M, Tassios PT, Daikos GL. Carbapenemases in *Klebsiella pneumoniae* and other Enterobacteriaceae: an evolving crisis of global dimensions. *Clin Microbiol Rev.* 2012;25(4):682-707.

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