Screening for bærertilstand af Vancomycin Resistente Enterokokker (VRE)

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5. Kuch A, Stefaniuk E, Ozorowski T, Hryniewicz W. New selective and differential chromogenic agar medium, chromID VRE, for screening vancomycin-resistant *Enterococcus* species. *J Microbiol Methods* 2009;77(1):124-126.
6. Peltroche-Llacsahuanga H, Top J, Weber-Heynemann J, Lutticken R, Haase G. Comparison of two chromogenic media for selective isolation of vancomycin-resistant enterococci from stool specimens. *J Clin Microbiol* 2009;47(12):4113-4116.
7. Delmas J, Robin F, Schweitzer C, Lesens O, Bonnet R. Evaluation of a new chromogenic medium, ChromID VRE, for detection of vancomycin-resistant *Enterococci* in stool samples and rectal swabs. *J Clin Microbiol* 2007;45(8):2731-2733.
8. Adler H, Oezcan S, Frei R. Vancomycin-resistant *Enterococci* of vanB genotype may pose problems for screening with highly selective media. *J Clin Microbiol* 2010;48(6):2323.
9. Mak A, Miller MA, Chong G, Monczak Y. Comparison of PCR and culture for screening of vancomycin-resistant *Enterococci*: highly disparate results for vanA and vanB. *J Clin Microbiol* 2009;47(12):4136-4137.
10. Marner ES, Wolk DM, Carr J et al. Diagnostic accuracy of the Cepheid GeneXpert vanA/vanB assay ver. 1.0 to detect the vanA and vanB vancomycin resistance genes in *Enterococcus* from perianal specimens. *Diagn Microbiol Infect Dis* 2011;69(4):382-389.
11. Zabicka D, Strzelecki J, Wozniak A et al. Efficiency of the Cepheid Xpert vanA/vanB assay for screening of colonization with vancomycin-resistant enterococci during hospital outbreak. *Antonie Van Leeuwenhoek* 2012;101(3):671-675.
12. Werner G, Serr A, Schutt S et al. Comparison of direct cultivation on a selective solid medium, polymerase chain reaction from an enrichment broth, and the BD GeneOhm VanR Assay for identification of vancomycin-resistant enterococci in screening specimens. *Diagn Microbiol Infect Dis* 2011;70(4):512-521.
13. Fang H, Nord CE, Ullberg M. Screening for vancomycin-resistant enterococci: results of a survey in Stockholm. *APMIS* 2010;118(5):413-417.

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