



Laser Capture Microdissection

DTU Veterinærinstituttet

The gastro intestinal microbiota

- The GI tract contains approx. 10^{14} bacteria
- 1 g of intestinal content can inhabit 10^{11} bacteria
- The GI tract contains more than 500 different bacterial species.
- Most of the GI bacterial species are anaerobic and not yet culturable.

DTU Veterinærinstituttet



P.A.L.M.

DTU
DTU Veterinærinstituttet

Bacterial identification

Abortion in sheep → Lung sample with lesions → In situ hybridization

16S rRNA gene = *Fusobacterium necrophorum*

Journal of Clinical Microbiology 2006 (44) 4537-4540
 Molecular and Cellular Probes 2006 (20) 330-336

DTU
DTU Veterinærinstituttet

DNA extraction by LCM

A

4200 μm^2

B

324 bp to 670 bp in size from formalin fixated samples

BioTechniques 2006 (39) 864-868

DTU
DTU Veterinærinstituttet

Bacterial interaction (gene expression)




in vivo analyses of *Actinobacillus pleuropneumoniae* expression during an infection

Snap-frozen tissue section
 ↓
 Immunohistochemistry
 ↓
 LCM
 ↓
 qRT-PCR


Journal of Microbiological Methods 2007 (69) 414-416

DTU
DTU Veterinærinstituttet

Bacterial-host interactions

- Transcriptome analysis on epithelial cell layer from gnotobiotic or conventional preterm piglets
- mRNA in situ hybridization on host cell and isolation of potential pathogenic or probiotic bacteria

 DTU Veterinærinstituttet

Laser capture microdissection and microbial ecology


The DOs and the limitations on using microdissection

DOs

- fast and reliable harvesting of thousands of cells, "single cells" and subcellular structures from histological sections, cytopins etc.
- selectively capture single neurons, fetal cells, sperm cells and chromosomes
- no limitations for your **genetic and proteomic analysis**
- **DNA or RNA analysis and microarray technology** with non contact pure samples
- preparation from paraffin and frozen sections, also directly from glass slides
- Isolate and cultivate **living cells** from fresh tissue

Limitations

- Its not a FACS or a microtweezer
- You need bacterial microcolonies to get enough material (RNA and DNA)
- LCM only makes sense if you have a solid or embedded material with substructures
- **Most known applications are for tissue sections and cell cultures!**

 DTU Veterinærinstituttet